

The Aaron Smoke Shark is available with two complete systems specifically designed for the OB/GYN:



The A950™ is the first High Frequency Desiccator with Cut. It features Cut, Blend, and Coagulation (60 watts), Fulguration/Desiccation (35 watts), and Bipolar (30 watts) modes with six presets. The easy-touch color-coded membrane pads, on the front panel, allow for quick selection of each setting. Power output may be increased or decreased via the 3-button handpiece or the rotary power-control.

The Aaron 950-G includes ...

- **A950 ElectroSurgical Generator**
 - 5 Extra Disposable Electrodes
 - 50 sharp and 50 blunt disposable dermal tips
 - A901, reusable single 3-button pencil
 - A910-P, 10-handpiece drapes
- **SEGYN Smoke Shark system**
 - SE01, Smoke Shark with pneumatic footswitch
 - SF18, 18-hour filter
 - 786T, 7/8" tube 6' long
- **A1202C, reusable grounding cord**
- **A1202, 5-solid grounding pads**
- **SERF, reducer fitting (attaches to vaginal speculum)**
- **A812-C, electroSurgical mobile stand**
 - A812, mobile stand
 - A812-E, height extender
 - A808-T, instrument tray
 - A812-BT, bottom tray
- **ESLK, electroSurgical loop kit**
 - ES09, 10 mm x 10 mm loop
 - ES11, 15 mm x 10 mm loop
 - ES12, 20 mm x 8 mm loop
 - ES13, 20 mm x 15 mm loop
 - ES16, 10 mm x 8 mm square



The A1250™ is a multipurpose electroSurgical generator incorporating automatic safety features like self-test circuits, audible tones, discreet outputs, isolated circuitry, and Bovie NEM™ (neutral electrode monitoring) pad sensing. It features 120 watts of cutting power, five modes of operation, and digital error detection to satisfy all surgical demands with safety, reliability, and convenience.

The Aaron 1250-G includes ...

- **A1250 ElectroSurgical Generator**
 - 5 Extra Disposable Electrodes
 - 2 Disposable Pencils
 - A1252C, Reusable Grounding Cord
 - ESRE-1, 5 Split Grounding Pads
- **SEGYN, Smoke Shark System**
 - SE01, Smoke Shark
 - SF18, Eighteen Hour Filter
 - 786T, 7/8" Tube 24" Long
 - FSSEP, Pneumatic Footswitch
 - SERF, Reducer Fitting (Attaches to Vaginal Speculum)
- **ESLK, ElectroSurgical Loop Kit**
 - ES09, 10mm x 10mm Loop
 - ES11, 15mm x 10mm Loop
 - ES12, 20mm x 8mm Loop
 - ES13, 20mm x 15mm Loop
 - ES16, 10mm x 8mm Square
- **ESMS-C, ElectroSurgical Mobile Stand with Smoke Evacuator Tray and Instrument Drawer**

Distributed By:



7100 30th Aveue North • St. Petersburg, FL 33710-2902
 U.S. Phone 1-800-537-2790 • Fax 1-800-323-1640
 Int'l Phone +1-727-384-2323 • Fax +1-727-347-9144
 www.boviemed.com • sales@boviemed.com

smoke shark™

Smoke Evacuator

A Sleek, Powerful, Smoke-Eating Machine!

- Easy To Operate
- Quiet and Effective Smoke Evacuation
- Long-life 18-hour Filter
- Low Cost Per Procedure
- Excellent Particle Capture
- Quality Made in the U.S.A.



The Aaron® Smoke Shark™ is designed to be easy to use and extremely effective at smoke and particle removal at the operative sight. While designed for use with all electroSurgical generators, it is also excellent for use in controlling laser smoke. The Smoke Shark's 18-hour filter combines four-stage filtration into one casing (Pre-filter, ULPA, Carbon, Post-filter), keeping cost per procedure down and minimizing waste. The four-stage filter captures particles down to 0.01 mm with 99.99995% efficiency giving you the confidence you need for a

safe work environment. The 24 ounces of carbon contained in the Smoke Shark filter, over twenty times the amount of a competing brand, does a superior job of odor elimination.

A variety of accessories are available to customize the Smoke Shark to your specific application. These include laser-resistant tubing, tubing reducers for attachment to smoke evacuating pencils or vaginal speculum, many other tubing configurations (sterile or non-sterile), a foot control, and a remote switch activator.

The Aaron Smoke Shark™ by Bovie®, is designed to be easy to use and extremely effective removing smoke and particles at the operative sight. While designed to work efficiently with all electrosurgical generators, it is also excellent for use in controlling laser smoke.

Style # SE01 / SE01-220

STANDARD EQUIPMENT unit with pneumatic footswitch, power cord, user's guide

DIMENSIONS

Height: 29.9 cm (11.8")
 Width: 17.8 cm (7.0")
 Depth: 46.4 cm (18.3")
 Weight: < 6.5 kg (< 14 lbs.)
 w/o filter

LINE VOLTAGE

SE01: 120 VAC
 SE01-220: 240 VAC

LINE FREQUENCY

50 – 60 Hz

WARRANTY

two years



TYPICAL SETUPS

ELECTROSURGERY

- SE01, Smoke Shark unit with pneumatic footswitch
- SF18, long-life 18-hour filter
- 786T 7/8" x 6' tube (case of 24)

or

- SETW 7/8" x 6' tube with wand (case of 24)
- FSSEP, Pneumatic footswitch
- SEAS, arm stand

OB/GYN

- SE01, Smoke Shark unit with Pneumatic footswitch
- SF18, long-life 18-hour filter
- 786T, 7/8" x 6' tube (case of 24)

or

- SETW, 7/8" x 6' tube w/wand (case of 24)
- FSSEP, Pneumatic footswitch

LASER

- SE01, Smoke Shark unit with Pneumatic footswitch
- SF18, long-life 18-hour filter
- 786T, 7/8" x 6' tube (case of 24)
- SELW, 7/8" x 8" laser-resistant wand (pack of 3)
- FSSEP, Pneumatic footswitch
- SEAS arm stand

SMOKE SHARK™ ACCESSORIES



SF18
 long-life
 18-hr. filter
 1/box



FSSE
 footswitch -
 will fit older
 smoke
 evacuators
 models -1/box



FSSEP
 Pneumatic
 footswitch
 included with
 the SE01
 1/box



SELKS
 laparoscopic
 tubing kit,
 sterile
 12/box



786T
 7/8" x 6' tube,
 non-sterile
 24/box



SEAS
 arm stand
 1/box



SERS
 remote switch
 activator
 1/box



SELVS
 laparoscopic
 tubing kit with
 valve, sterile
 12/box



786TS
 7/8" x 6' tube,
 sterile
 24/box



SETW
 7/8" x 6' tube
 with wand &
 tip, non-sterile
 24/box



SEFC
 fluid collection
 canister 1200cc
 6/box



SERF
 reducer fitting
 7/8" - 1/4" x 24"
 10/box
 [fits optional
 vaginal speculum]



SELW
 laser-resistant
 wand 7/8" x 8",
 non-sterile
 3/box



SETWS
 7/8" x 6' tube
 with wand &
 tip, sterile
 24/box



SERH
 canister ring
 holder - 1/box



SERFS
 reducer fitting
 7/8" - 1/4" x 24"
 sterile - 10/box
 [fits optional
 vaginal speculum]



SELWS
 laser-resistant
 wand 7/8" x 8",
 sterile
 10/box



SELM
 laser mask
 50/box



SEWP
 wall plate to
 hold canister
 1/box



SEVL
 7/8" x 6' tubing
 reduced to 3/8"
 for ValleyLab &
 IEC pencil,
 sterile - 5/box