

OPTIMAL CONTROL BETTER PROTECTION

THE HEALTH HAZARDS OF ELECTROSURGERY
IN THE PHYSICIAN'S OFFICE



THE Bovie® SMOKE EVACUATION SOLUTION

THE SMOKE PLUME CHALLENGE


80% of surgeries worldwide use electrosurgery
Smoke plume is created in 95% of these procedures

biological, physical and chemical hazards

limitations



Multiple commissions and associations have established guidelines for addressing surgical smoke, with some states having passed legislation to make smoke evacuation mandatory



THE BOVIE® SMOKE EVACUATION SOLUTION




INTRODUCING THE NEW BOVIE® SMOKE SHARK® III

- Uniquely designed for the physician's office
- Specifically intended for shorter (<15 minutes) procedures
- Reduces biological, physical and chemical hazards present in electrosurgical smoke plume
- Minimal noise level with efficient filtration
- Filter life tracking technology
- Hand-operated activation for maximizing filter life

4-Stage ULPA Filtration

99.999%

efficient in trapping 0.1-0.2 μm particulates and micro-organisms



Bovie® Smoke Shark® III & High Frequency Dessicator (shown to scale)



Smoke Protection

- 3/8" port filtration tubing connection
- 2-year warranty

Size

- Lightweight and compact for simplified operation, storage and portability; 9.4" (23.9 cm) wide x 7.1" (18 cm) tall x 6.1" (15.5 cm) deep
- Ultra-quiet motor (<55 dB) provides powerful vacuum for efficient smoke evacuation

Economical

- Accommodates SharkSkin smoke evacuation pencil adapter to work with most pencils
- Multi-procedure filter (SF06) provides up to 6 hours of performance
- Smart Shark remote activation switch maximizes filter life

SMOKE SHARK® III

UNIQUELY DESIGNED FOR THE PHYSICIAN'S OFFICE



Feature	Smoke Shark III	Smoke Shark II
Smaller footprint for ease of use, transportation and storage in facilities where space is at a premium	Y	
Filter life designed for multiple, shorter procedures	Y	
Lower cost for facilities with limited financial resources	Y	
ESU Sensor Activation switch included to maximize filter life	Y	
Ultra-quiet vacuum to gather smoke plume (<55dB)	Y	Y
ULPA filter to remove and trap 0.1-0.2µm particles and micro-organisms at 99.999% efficiency	Y	Y

SYSTEM INCLUDES

Multi-procedure 6-hour filter	Smart Shark remote activation switch for filter life optimization	SharkSkin smoke evacuator pencil adapter (with 10' tubing)
		
Included (2-Pack Sold Separately)	Included	Individual unit included (Box of 20 Sold Separately)

Bovie® Smoke Shark® III System

CODE	DESCRIPTION
SE03-120	Smoke Shark® III 120V unit
SE03-220	Smoke Shark® III 220V unit

Accessories

CODE	DESCRIPTION
SF06	Smoke Shark® III Filter 2-Pack (non-sterile)
SERS2	Smart Shark remote activation switch
SS95	SharkSkin smoke evacuation pencil adapter with 3/8" x 10' Tubing (20/bx, packaged sterile individually)

CALL YOUR SYMMETRY SURGICAL SALES REPRESENTATIVE, LEARN MORE AT SYMMETRYSURGICAL.COM OR CALL 1-800-251-3000.

3034 Owen Drive | Antioch, TN 37013 | 615.883.9090 | 800.251.3000 | customerservice@symmetrysurgical.com | symmetrysurgical.com

©2020 Symmetry Surgical, Inc. All rights reserved. SM0093
08.07.20V1