

Midmark 550 Stretcher - Refurbished



Stock Photo – Options May Vary

The Midmark 550 Stretcher is a Versatile, Full-Functioning Stretcher that is perfect for the Operating Room, Emergency Room, Post-Operative Care Unit or Transport. In addition, the Midmark 550 has a Pneumatic Assist Fowler (backrest) 90°, Hydraulic Height Adjustment (foot pedal on foot end of stretcher) 24" to 33", and Trend/Reverse Trend.

The Midmark 550 has a Weight Capacity of 500 lbs., a Central Brake, Four Wheel Locking and Steering, 8" Casters for easy maneuverability and control, along with Fold-Down side Rails. Midmark Stretcher 550 comes with O2 Tank Holder, New Pad, and 6 Months Parts Warranty. See below for more information.

The Midmark 550 Stretcher comes completely refurbished with a comprehensive warranty and a new pad.

- Completely Refurbished
- Trendelenburg/Rev Trendelenburg
- Head Up/Down (Pneumatic Assist)
- Height Up/Down (Hydraulic)
- Fold Away Side Rails
- 4-Wheel Locking & Steering
- O2 Tank Holder
- 500 lbs Weight Capacity
- New Pad
- 6 Months Parts Warranty

Specifications:

- Overall Dimensions: 29.5"W x 83"L
- Height: 24" - 33"H
- Height Up/Down (Hydraulic)
- Side Rails: Fold Away
- 4-Wheel Locking & Steering
- 500 lb Weight Capacity
- 6 Months Parts Warranty

***Brochure and images are a general overview of the product. Actual product options may vary.*