



www.venturemedical.com

Venture Medical ReEquip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Steris/Amsco Quantum Surgical Lights - Refurbished



Stock Photo – Options May Vary

The Amsco Quantum surgical lighting system is designed to provide exceptional illumination in almost any surgical setting. The 24" light head is easily maneuvered, and the auto re-lamping feature means that you should never be without light, even if a bulb fails during surgery.

The large, pre-focused beam of light is pre-focused, and yields an incredible 10,000-foot candles (107,600 lux). Impressive natural color rendition at 4200 Kelvin and a color rendering index (CRI) of 92 provides are also key features of this system.

The Steris Amsco Quantum surgical light can come in a single or dual head configuration. Also included are the wall controller and transformer. The surgical lights come fully refurbished and includes a warranty.

- Flexible, Versatile Lighting
- Excellent Depth of Illumination (26")
- Automatic Re-Lamping
- Color temperature 4200k
- Includes 360 Degree Swivel Arm Support
- Intensity: 10,000 ft candles (107,600 lux)
- Wall Mounted Intensity Controller
- Color Rendering Index: 92
- Focusable Light Head
- Excellent Shadow Control
- Quartz Halogen Bulbs
- Refurbished

Specifications:

- Halogen Bulb: 22 V / 220 W
- Intensity Pattern (large): 7.5"
- Intensity Pattern (small): 6.75"
- Clearance (rotation): 14.5'
- Intensity: Up to 10,000 Ft Candles
- Diameter of Light Head: 24"
- Power Requirement (volts/Hz):
 - 110-120V, 50/60 Hz
 - 230-240V, 50/60 Hz
- Power Requirements (amps):
 - 7A (dual head)
 - 3.5A (single head)

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReEquip, Inc.