



## Preferred Care® Bed-Recliner

### SPECIFICATIONS & ORDERING INFORMATION

**AR981642** Lumex® Preferred Care® Bed-Recliner; includes:

1. Sleep Surface and Overlay
2. Lit Patient/Caregiver Hand Pendant
3. One-Step Locking System
4. Wall Saver
5. Battery Back-up
6. Under Bed Light
7. Edema Ratchet

**AROL100** Lumex Overlay

**ZL872067** Lumex® Preferred Care® Bed-Recliner Low Profile Head & Foot Boards

**ZA84200** Half Length Assist Device Warm to the Touch

**ZA85400** Pivot Assist Bar

**650-1001-000** Assist Bar for Lumex® Preferred Care® Bed-Recliner

### Specifications & Dimensions

- Mattress Deck Length: 80"
- Mattress Deck Width: 35"
- Overall Bed Length Without Wall Saver: 85.5" w/Head & Foot Boards
- Overall Bed Width: 41.0" w/Assist Device and/or Assist Bar
- Maximum Deck Height: 30"
- Minimum Deck Height: 7"
- Maximum Deck Height w/Sleep Surface & Overlay: 37.5"
- Minimum Deck Height w/Sleep Surface & Overlay: 14.5"
- Safe Workload (including resident, support surface & accessories): 600 lbs  
Max Weight capacity evenly distributed.

The Lumex® Preferred Care® Bed-Recliner is built to conform to UL 60601-1 and IEC601-2-38; CAN/CSA C22.2 No. 601.1. Multiple Patents Pending

### Warranty Information

- Bed Frame: Limited Lifetime Warranty
- Welds: Limited 15 Year Warranty
- DC Electronic Components: Limited 3 Year Warranty
- Head/Foot Boards: Limited 1 Year Warranty
- Sleep Surface & Overlay: Limited 2 Year Warranty

The warranties contained herein contain all the representations and warranties with respect to the subject matter of this document, and supersede all prior negotiations, agreements and understandings with respect thereto. The recipient of this document hereby acknowledges and represents that it has not relied on any representation, assertion, guarantee, warranty, collateral contract or other assurance, except those set out in this document.

Please see the complete warranty statement at [www.grahamfield.com](http://www.grahamfield.com)



GF Health Products, Inc. Atlanta, Georgia 30360  
800.347.5678 [www.grahamfield.com](http://www.grahamfield.com)



Specifications are subject to change without notice.  
© 2012, Lumex and SilverSolutions are registered trademarks of GF Health Products, Inc. All rights reserved. GF1000051RevB12

### Patient Comfort Combined With Ease of Maintenance and Durability

The Sleep Surface and Lumex Overlay has been designed to maximize the comfort of the patient in all positions.

Four corner straps secure the Lumex Overlay to the Sleep Surface. The Sleep Surface is constructed of durable foam covered in Silvertex Fabric featuring Silverguard® Silver Ion Technology.

### Sleep Surface Construction:

3.5" High Density Base Foam  
1.5" Top Wrap Comfort Layer Foam  
Upholstery and Trim: Silvertex Marine Blue: Meets CA TB 133

### Lumex Overlay Construction:

1.5" High Density Base Foam  
Top Foam: Visco-Elastic 1.5" Foam  
Top Cover: Fluid-Resistant, Five Way Stretch Fabric  
Bottom Cover: Fluid-Resistant, Non-Skid Material  
Elastic Attachment Straps at Each Corner

For maximum durability and infection control, the epoxy coated frame features **SilverSolutions**. Integrated directly into the finish, **SilverSolutions** kills more than 99.9% of odor and stain causing bacteria. **SilverSolutions** is non-allergenic, safe for the environment and lasts for the expected life of the bed.



## Preferred Care® Bed-Recliner

COMFORT • VERSATILITY • SAFETY



### Designed For:

- Dialysis
- Nocturnal Dialysis
- Oncology
- Cardiac Care

- Same Day Surgery
- Post-Operative Recovery
- Any Time-Intensive Procedure





# Preferred Care® Bed-Recliner

COMFORT • VERSATILITY • SAFETY

Designed for Dialysis, Nocturnal Dialysis, Oncology and any time intensive health-care procedures, the Lumex® Preferred Care® Bed-Recliner combines maximum patient comfort with ease of use by the caregiver. Unlike recliners that are converted into a semi-flat sleeping surface, the Lumex® Preferred Care® Bed-Recliner offers the comfort of a true bed/sleep surface with an almost infinite number of seated/reclined positions. A specially designed Sleep Surface and Flexible Overlay provides proper support and extreme comfort in any position.

Ultra Quiet DC Motors allow the Preferred Care® Bed-Recliner to quickly move from a low position of 7.0" Deck Height to 30" Deck Height (without Sleep Surface and Overlay). This range of motion allows for easier care by the caregiver in the full up position and helps prevent patient injuries in the low position. The bed deck can be lowered to a height of 7.0" from the floor. Timesaver Technology-fast motors save time when changing positions/ height of Bed-Recliner.



Shown with Sleep Surface and Overlay

### STANDARD FEATURES:

- Full and Reverse Trendelenburg Positions
- Auto CPR Function
- Infinite Positioning of Head and Foot Sections
- Adjustable Leg-Lift System to Improve Circulation
- No Pinch Points Within 200mm of Bed Perimeter
- Rolling Clearance of 20mm Over Thresholds
- Splash-Proof Electronics
- Wall Bumper
- Edema Foot Ratchet
- Under Bed Light
- Battery Back Up

### Patient Comfort & Caregiver Features/Options for the Lumex® Preferred Care® Bed-Recliner



The Preferred Care Bed-Recliner features a Patient Hand Pendant (Backlit) and a Low-Profile Footboard with an Embedded Caregiver Control Panel. The Embedded Caregiver Control Panel can be used to "lockout" the Patient Hand Pendant preventing unauthorized use of the features of the Lumex® Preferred Care® Bed-Recliner.

The safety inspired **embedded staff control** in footboard is:

- Angled for better visibility
- Features Two-stage lockouts – lockout of pendant use and lockout of embedded staff control
- Quick lockout buttons lock multiple functions simultaneously
- SafeMode™ automatically locks out high/low feature (unless disabled), Patents Pending

**Standard Under Bed Light** for night-time safety.

#### Edema Foot Ratchet

Provides quick positioning for improved circulation in the patient's legs. Standard on the Lumex® Preferred Care® Bed-Recliner.



#### Cardiac Position

The Cardiac Position is primarily designed for residents with heart and respiratory problems to whom a semi-erect position brings great relief.

The Cardiac Position helps to elevate the patient from the waist up, which provides relief to the lungs, increases circulation, and assists the patient in coming from a fully prone position to a sitting position without causing harm or undue strain.

#### Two-Pedal Locking System

Standard on the Lumex® Preferred Care® Bed-Recliner, the clearly indicated Red and Green Coded Locking System allows the bed to be locked at any height. Easily accessible from the foot of the Preferred Care® Bed-Recliner, the bed can be locked or unlocked with just one step on the red coded locking pedal. To unlock the bed, simply step down on the green coded locking pedal.

#### Battery Assist

In case of a power failure, the Preferred Care® Bed-Recliner motors include a battery assist to power the basic functions of the bed.



**Head-Up Position:** Ideal for Resting and Treatment Initiation



**Sleep Position:** Ideal for Sleeping and Patient Care



**Recline/TV Position:** Ideal for Resting, TV & Activities



*The almost infinite positioning and height adjustment capability of the Preferred Care® Bed-Recline makes it ideal for patient transfer from a wheelchair, stretcher or another bed. The ability to have equipment at the same height helps to make patient transfers safer for both the patient and caregiver.*

**Optional Half Length Assist Device and Assist Bars** provide a sturdy and secure hand-hold to assist residents in and out of bed. Assist devices are removable without tools and meet FDA Entrapment Guidelines.



#### Half Length Assist Device-Warm to the Touch (ZA84200)

- Warm-to-the-touch half-length assist device and pivot assist bars:
- Move out of the way for lateral transfers
- Pivot away from egress area
- Provide toe clearance even in lowest position to prevent staff injury
- Hand pendant mounts facing inward or outward and away from hand grab area



#### Assist Bar (650-1001-000)

Can be mounted without tools at the head section of the Lumex® Preferred Care® Bed-Recliner or the foot section of the bed based on the patient's requirements. Provides a sturdy hand-hold to assist patients in and out of bed.