

# Elite Comfort Recliner 5900



A gem of a chair for its superior comfort and durability, crafted with style

Attention to detail, concentration on comfort and quality, delivered at the best value; standards that Winco recliners are known for. Winco's Elite Comfort Recliner rises above the competition and are among the best-selling chairs on the market today.



**HEAT  
&  
MASSAGE**  
(Optional)

## KEY FEATURES

- Black 3" locking casters
- Built-in headrest
- Fully upholstered arm side panels

Casters Shown are optional,  
Standard Casters are black.

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

[www.wincomfg.com](http://www.wincomfg.com)



# PRODUCT SPECIFICATIONS

## ELITE COMFORT RECLINER 5900

- 1 Push handle on back for easy maneuverability.
- 2 Built-in headrest.
- 3 Wide elastic bands on back frame provide enhanced comfort and support.
- 4 Non-porous, fold-down side tray contains small spills. STANDARD on right side.
- 5 Black 3" Nylon locking casters provide ease of maneuverability (White casters shown, optional)










### 5900

Seat Height . . . . .	18.5" (47cm)
Seat Width . . . . .	21" (53.3cm)
Seat Depth . . . . .	20" (50.8cm)
Overall Height . . . . .	46" (116.8cm)
Overall Width* . . . . .	33.5" (85.1cm)*
Back Height Above Seat . . . . .	32.5" (82.6cm)
Overall Length (Upright) . . . . .	41" (104.1cm)
Overall Length (Reclined) . . . . .	63" (160cm)
Arm Height . . . . .	26" (66cm)
Weight . . . . .	97 lbs. (44kg)
Weight Capacity . . . . .	275 lbs. (124.7kg)

\*Add 9" to chair width for each raised side table

SPECIFICATIONS

OPTIONS

 IB-133 <b>TB</b>	 Heat <b>HT</b> Heat & Massage <b>HM</b>	 Armrest Covers <b>AC00</b>	 Side Cushion <b>SC00</b>
 IV Pole <b>IV</b> (Left rear only)	 Leg Extensions <b>SE00</b>	 Left Side Tray <b>LT</b> (additional)	

### STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

COLORS

### OPTIONAL COLOR UPGRADES

(ADDITIONAL CHARGES APPLY)



### SILVERTEX VINYL

(ADDITIONAL CHARGES APPLY)

