

Adjustable Back Treatment Table

8570



ISO 9001-2008 Certified

Superior durability, outstanding features, and an affordable price

Our adjustable back tables offer the durability and range of standard features you expect at a price that you'll love. One look, and you'll see the difference.



8570

KEY FEATURES

- Gas spring adjustable backrest
- Durable foam construction
- Sturdy steel leg frames
- Paper holder
- Adjustable leveling feet



8570-SH



8570-C1



8570-D1

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

ADJUSTABLE BACK TREATMENT TABLE 8570

- 1 Adjustable headrest to enhance patient comfort. Angle range: 6-57 degrees
- 2 One-touch handle easily raises headrest to desired angle.
- 3 Standard paper holder.
- 4 High strength powder-coated steel frame.
- 5 High impact leveling feet ensure a sturdy table on almost any hard surface.
- 6 Heavy gauge expanded vinyl & medium density foam for superior comfort.



SPECIFICATIONS

8570

Overall Height*	30" (76cm) STANDARD*
Overall Width	28" (71cm)
Overall Length	72" (183cm)
Adjustable Backrest Range	6 - 57 degrees
Weight**	91 lbs. (41.3kg)**
Weight Capacity	400 lbs. (181.4kg)

*Special height legs available in 25" or 36"
 **Weight varies depending on table options ordered.
 Add 25 lbs. for 8570-SH

OPTIONS



TB-133 TB



*Special Height Legs



Shelf SH



Cabinet C1



Drawer D1



Matching Pillow MPOO

STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



COLORS