

## Overbed Table, Plastic Top

HCPCS:



- Plastic top is durable, easy to clean and has a raised edge to contain spills
- Table surface won't crack or peel and includes two areas for placing drinks
- Tabletop can be raised or lowered in settings between 29" - 42" with light upward pressure
- Tabletop locks securely when height adjustment handle is released
- Swivel casters allow for easy maneuverability
- Strong "H" base provides security and stability
- Comes in only four parts for quick and easy assembly

**PRODUCTS****Description**

| PRODUCTS | Description                             |
|----------|---|
| 13069    | OVERBED TBL W/PLASTIC TOP SLV VEIN, 1/E |

**SPECIFICATIONS****Description**

| SPECIFICATIONS             | Description       |
|----------------------------|-------------------|
| Base Dimensions            | 29" (L) x 13" (W) |
| Carton Shipping Weight     | 22.93 lbs         |
| Height Adjustment          | 29" - 42"         |
| Tabletop Surface Dimension | 28" (L) x 15" (W) |
| Warranty                   | Limited Lifetime  |
| Weight Capacity            | 50 lbs            |