

## Overbed Table, Plastic Top

HCPCS:



- Plastic top is durable, easy to clean and has a raised edge to contain spills
- Table surface won't crack or peel and includes two areas for placing drinks
- Tabletop can be raised or lowered in settings between 29" - 42" with light upward pressure
- Tabletop locks securely when height adjustment handle is released
- Swivel casters allow for easy maneuverability
- Strong "H" base provides security and stability
- Comes in only four parts for quick and easy assembly

PRODUCTS	Description
13069	OVERBED TBL W/PLASTIC TOP SLV VEIN, 1/E
SPECIFICATIONS	Description
Base Dimensions	29" (L) x 13" (W)
Carton Shipping Weight	22.93 lbs
Height Adjustment	29" - 42"
Tabletop Surface Dimension	28" (L) x 15" (W)
Warranty	Limited Lifetime
Weight Capacity	50 lbs