

Amsco M70 Warming Cabinet - Refurbished



Stock Photo – Options May Vary

The Amsco M70 Blanket Warmer is a cost-effective solution for blanket and fluid warming needs in your medical facility.

The M70 features analog controls for the top compartment to adjust temperature; the bottom cabinet maintains a constant 120° for blankets. This unit has both visual and audible alarms in case of over-temperature detection.

The Amsco M70 warming cabinet comes completely refurbished with a comprehensive warranty and 1 shelf.

- Dual Compartment
- Cost Effective
- Top and Bottom Heating
- Adjustable Temp in Upper Compartment (Dial)
- 120° Temp in Bottom Compartment
- Standard 120V Plug
- Stainless Steel Construction
- Over Temp Alarms (Audible and Visual)
- Includes Blanket Rack and 1 Bottom Compartment Shelf
- Completely Refurbished

Specifications:

- Dims: 26.5" D x 30" W x 74.75" H
- Power: 120V 7.5 Amps
- Standard AC Wall Plug
- Approx Weight: 320 lbs
- Temp Range: 95° - 150°
- Temp Range Upper: 95°-150°
- Temp Lower (Constant) – 120°

***Brochure and images are a general overview of the product.
Actual product options may vary.*