

## Baxter FloGard 6201 Infusion Pump - Refurbished



Stock Photo – Options May Vary

The Baxter FloGard 6201 single channel all-purpose volumetric infusion pump uses standard IV administration sets. The 6201 features automatic restart after an occlusion clears and IV set free-flow protection.

The FloGard 6201 infusion pump can deliver a wide variety of fluids over a broad range of infusion rates. The 6201 operates on a rechargeable battery that lasts approx. 6 hours, recharge takes approx. 8 hours.

The Baxter FloGard 6201 comes completely refurbished with a comprehensive warranty.

- Battery Operation
- Flow Check On/Off Display
- Slide Clamp Option
- 3 Levels of Sensitivity
- Audible Alarm
- Dose Checking
- Electric Memory
- Completely Refurbished

### Specifications:

- Overall Dims: 7.9" x 5.1" x 11.4"
- Battery: 6 hrs at 1,400 mL/hr
- Battery Recharge: 80% within 8 hrs.
- Battery: 12 V, Lead Acid
- Communication Port: R5232
- Weight: 11.6 lbs
- Power Requirements:
  - 115 VAC, 60Hz, 30W

*\*\*Brochure and images are a general overview of the product. Actual product options may vary.*