

Baxter FloGard 6301 Infusion Pump - Refurbished



Stock Photo – Options May Vary

- Battery Operation
- Flow Check On/Off Display
- Slide Clamp Option
- 3 Levels of Sensitivity
- Audible Alarm
- Dose Checking
- Electric Memory
- Completely Refurbished

The Baxter FloGard 6201 dual channel volumetric infusion pump uses standard IV administration sets. The 6301 features automatic restart after an occlusion clears and IV set free-flow protection. The two separate pump channels allow it to do the work of two conventional pumps, resulting in space savings.

The FloGard 6301 infusion pump can deliver a wide variety of fluids over a broad range of infusion rates. The 6301 operates on a rechargeable battery that lasts approx. 6 hours, recharge takes approx. 8 hours.

The Baxter FloGard 6201 comes completely refurbished with a comprehensive warranty.

Specifications:

- Overall Dims: 13" x 5.1" x 11.4"
- Battery: 6 hrs at 1,400 mL/hr
- Battery Recharge: 80% within 8 hrs.
- Battery: 12 V, Lead Acid
- Communication Port: R5232
- Weight: 17.9 lbs
- Power Requirements:
 - 110/120 V, 60Hz

***Brochure and images are a general overview of the product.
Actual product options may vary.*