

## Baxter AS50 Syringe Pump - Refurbished



Stock Photo – Options May Vary

Syringe Types:	B-D: 1, 3, 5, 10, 20, 30, 60 mL plastic Monoject: 1, 3, 6, 12, 20, 35, 60 mL Terumo: 1, 3, 5, 10, 20, 30 60 mL
Flow Rate Accuracy:	±3% of full-scale plunger travel
Volume Accuracy:	±3% For volume infusions – 1.3%, or .007” of travel
Flow Rate Range:	0.01 mL/hr to 438 mL/hr, depending on the syringe
Deliverable Volume:	Full syringe volume for 1-60 mL syringes
Data Display:	Self-prompting, multi-field LCD
Status Display:	Nine LED array
Battery Life:	5 hours of operation at 100 mL/hr 12 hours of operation at 2mL/hr using a 60 mL syrin following a charge of not less than 16 hours

The Baxter AS50 syringe pump is a compact unit that features a built-in, configurable drug library, syringe plunger detection, a backlit display, adjustable occlusion specialties and selectable modes of drug delivery.

The AS50 also features rubber end bumpers to help prevent damage from accidental drops. Additionally, the syringe pump can be mounted directly to a standard IV pole with the included IV pole mounting bracket.

Included with the unit: New Battery, IV Pole Clamp, Charger and Operators Manual

The Baxter AS50 comes completely refurbished with a comprehensive warranty.

### Specifications:

- Overall Dims: 3.4” W x 2.6” D x 10” H
- Weight: 3.2 lbs
- Power: AC: 105-125V (Battery Charger)
- Battery: Internal Nickel Cadmium

*\*\*Brochure and images are a general overview of the product.  
Actual product options may vary.*