

## Baxter InfusOR Syringe Pump - Refurbished



Stock Photo – Options May Vary

The Baxter Infus O.R. Pump is a syringe infusion device that will aid in the administration of many intravenous agents given during anesthetic procedures. In conjunction with the Smart Label System, it provides for the convenient delivery of narcotics, muscle relaxants, vasoactives and other drugs routinely given during anesthesia. The Smart Label System is a set of drug specific labels that attach to the front of the INFUS O.R. Pump. Each Smart Label converts the INFUS O.R. Pump into a specific delivery system for the indicated drug with appropriate rate selections and delivery units. By simply changing the specially coded Smart Label, the INFUS O.R. Pump is automatically changed to a specific delivery system for a particular intravenous agent.

Included with the unit: New Batteries, IV Pole Clamp, and Operators Manual

The Baxter InfusOR comes completely refurbished with a comprehensive warranty.

- User Friendly Design
- Modify dose without multiple reprogramming steps.
- Adjust settings easily at any time during infusion.
- Change dial settings with a single hand for an immediate response.
- Save calculation time.
- Minimize risk of calculation error.
- Change Infusion Parameters/Rates with Smart Labels
- Smart Labels are magnetically coded/drug specific.
- Select from the available SMART Label library.
- Ideal for administering drugs requiring multiple dosage adjustments.
- IV Pole Mountable

### Specifications:

- Dims: 9.2"H x 4.5"W x 2.0"D
- Weight: 2 lbs (Including Batteries)
- Accuracy:  $\pm 3\%$
- Syringes: Monoject or BD
- Battery Life: Approx 150 hours

### Smart Labels Available:

Propofol	Sufentanil
Dopamine	Mivacurium
General mL/hr	Cisatracurium
Fentanyl	Remifentanyl
Sufentanil	Vecuronium