

## Gendron 890 Stretcher - Refurbished



*Stock Photo – Options May Vary*

Our Completely Refurbished Gendron 890 Stretcher has all of the features you need at a great value! The Easy Access Gendron 890 Stretcher features Trendelenburg and Reverse Trendelenburg, Pneumatic Assist Head Up/Down, easy Height Up/Down, with conveniently located foot pedals on either side.

**Additional Gendron 890 Features:** The Central steering/braking system can be operated from either end/side. Fowler assisted Backrest raises 90° - when in the upright position, conveniently allows access to the storage area, which houses an IV Pole (sold separately), along with a security cable and an O2 Tank Well. Full-length Siderails fold down and under out of the way to allow easy access. The 890 has large 8" all Surface Carpet Casters - four-wheel central casters lock and swivel/steering lock.

- Trendelenburg/Rev Trendelenburg
- Head Up/Down (Pneumatic Assist)
- Height Up/Down
- Fold Under Siderails
- 4 Wheel Locking
- Steering
- O2 Tank Well
- New Pad
- 6 Months Parts Warranty

### Specifications:

- Fowler assisted Backrest: 0- 90°
- 28"x76" patient surface
- Caster Size: 8"
- 4 Wheel Locking Swivel Casters
- 6 Months Parts Warranty

*\*\*Brochure and images are a general overview of the product. Actual product options may vary.*