

Hill-Rom Transtar P-8000 Stretcher - Refurbished



Stock Photo – Options May Vary

The Hill-Rom Transtar Stretchers are perfect for emergency departments, transport, PACU, procedures, or anywhere else a versatile stretcher is needed. Higher weight capacity stretchers are also available.

Standard features include a lift assist backrest (Fowler), Trendelenburg, reverse Trendelenburg, and height adjustment. The 500 lbs weight capacity can handle large patients with ease.

The Hill-Rom Transtar stretcher comes completely refurbished with a comprehensive warranty and a new pad.

- Gas Spring Assist Backrest
- Collapsible Side Rails
- Central Brake and Steer
- \pm Height, \pm Trendelenburg
- Built-in O2 Tank Holder
- 500 lbs Weight Capacity
- Includes New 4" Pad
- Completely Refurbished

Specifications:

- Overall Length: 83"
- Overall Width: 26.5" or 30.25"
- Patient Surface Width: 26"
- Height Range: 21.5" – 35"
- Weight Capacity: 500 lbs
- Trend/Rev Trend: \pm 18°
- Pad Thickness: 4"
- Optional Knee Gatch

**Specifications are approximate.*

**Specifications will vary by model and features selected.*

***Brochure and images are a general overview of the product.
Actual product options may vary.*