

## Medfusion 3500 Syringe Pump - Refurbished



Stock Photo – Options May Vary

The Medfusion™ 3500 Syringe Pump is an intuitive system designed to enhance safety while maintaining simplicity. For more than 15 years, the Medfusion™ Syringe platform has been the technology leader in syringe pumps with its innovations in syringe recognition, drug libraries and dose limits. Support for up to 16 different care areas allows for transport of the Medfusion™ Syringe Pump into any department within the facility. As measured by the ISMP and other regulatory qualifications, the Medfusion™ Syringe Pump is truly a "smart" pump.

Included with the unit: New Battery, IV Pole Clamp, Power Cord and Operators Manual

The Medfusion 3500 comes completely refurbished with a comprehensive warranty.

- Easy to Use
- Programmable Dose Limits
- Rapid Occlusion Detection
- Post-Occlusion Bolus Reduction
- 13 Delivery Options
- Various Delivery Modes
- Connectivity with EMR
- Up to 10 hrs Battery Capacity
- Accepts Most Syringes up to 60cc
- IV Pole Mountable
- Completely Refurbished
- Warranty

### Specifications:

- Dims: 5.1"H x 10.5"W x 5.5"D
  - Weight: 3.8 lbs (pump only)
  - Accuracy: +/- 2%
  - Flow Rate: 0,01 ml/hr to 130 ml/hr
- \*\*Depends on Syringe\*\*

*\*\*Brochure and images are a general overview of the product.  
Actual product options may vary.*