



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Midmark 535 Stretcher - Refurbished



- Trendelenburg/Reverse Trendelenburg
- Back Up/Down (Pneumatic Assist)
- Height Up/Down (Foot Pump at Foot End of Stretcher)
- Fold-Under Side Rails
- 4 Wheel Locking
- Steering
- 350 lbs Weight Capacity
- New Pad
- 6 Months Parts Warranty

The Cost Effective, Completely Refurbished, Midmark 535 Stretcher features Include: Pneumatic Assist Fowler (backrest) 90°, Hydraulic Height Adjustment (foot pedal on foot end of stretcher), Trend/Reverse Trend, Central Brake and Fold-Down Side Rails. Stretcher has a 350 lbs Weight Capacity.

Additional Features: The 535 has 8" Casters for easy maneuverability and control, 12" x 63 ¾" Side Rails, and each Midmark Stretcher comes with a New Pad and Warranty. See below for more information.

The Midmark 535 Stretcher comes completely refurbished with a comprehensive warranty and a new pad.

Specifications:

- Overall Dimensions: 84" x 33"
- Back Up/Down (Pneumatic Assist); 0° - 90°
- Height Up/Down: 25¾" - 35¾"
- Fold-Under Side Rails
- Brakes: 4 Wheel Locking
- Casters: 8"
- Weight Capacity: 350 lbs.
- 6 Months Parts Warranty
- Weight: 268 lbs.

VENTUREMEDICAL.COM

© 2007-2025 Venture Medical ReQuip, Inc.