



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Midmark 545 Stretcher -Refurbished



- Trendelenburg/Rev Trendelenburg
- Head Up/Down (Pneumatic Assist)
- Height Up/Down (Hydraulic)
- Fold Away Side Rails
- 4 Wheel Locking
- Steering
- O2 Tank Holder
- 350 lbs Weight Capacity
- Includes New 4" Pad
- Completely Refurbished

The Midmark 545 Stretcher is a versatile, full-functioning stretcher that is perfect for the operating room, emergency room, post-operative care unit or transport. Features include pneumatic assist fowler, hydraulic height adjustment and trend/reverse trend.

Additional features include: four wheel locking and steering, 8" casters for easy maneuverability and control, a 350 lbs weight capacity along with fold-down side rails.

The Midmark 545 Stretcher comes completely refurbished with a comprehensive warranty and a new pad.

Specifications:

- Weight Capacity: 350 lbs.
- Overall Length: 83"
- Overall Width: 33.5"
- Height Range - Low: 24"
- Height Range - High: 33"
- Fowler Backrest: 0° - 90°
- Trendelenburg: +12°
- Reverse Trendelenburg: -7°
- Patient Surface: 77" x 31"
- Side Rails: 15" x 63 5/8"
- Casters: 8"
- Pad Thickness: 4"

VENTUREMEDICAL.COM

© 2007-2025 Venture Medical ReQuip, Inc.