



www.venturemedical.com

Venture Medical ReQuip, Inc.  
6008 Bonacker Drive  
Tampa, FL 33610  
800-627-3215 | 813-740-8640  
Fax: 800-466-1251  
sales@venturemedical.com

## Stryker 1067 Eye/Head/Neck Stretcher - Refurbished



\*\*Wrist Rest (pictured) not included.

- Collapsible Side Rails
- Central Brake and Steer
- $\pm$  Height,  $\pm$  Trendelenburg
- Crank Operated Backrest
- Built-in O2 Tank Holder
- 500 lbs Weight Capacity
- Includes New 2" Pad
- Completely Refurbished

The Stryker 1067 procedural stretcher features a tapered head rest for performing procedures on the eye, head and neck area.

Features include: hand crank backrest (fowler), trendelenburg, reverse trendelenburg, and height adjustment. The 500 lbs weight capacity can handle large patients with ease.

The Stryker 1067 stretcher comes completely refurbished with a comprehensive warranty and a new pad.

### Specifications:

- Overall Length: 84"
- Overall Width: 28.25"
- Patient Surface Width: 24"
- Height Range: 24.5" – 37.5"
- Weight Capacity: 500 lbs
- Fowler Angle: 0° - 90°
- Trend/Rev: +12° to -18°
- Wheel Diameter: 8"
- Pad Thickness - 2"

\*Please note that the head rest is tapered, but does not articulate.

**VENTUREMEDICAL.COM**

© 2007-2025 Venture Medical ReQuip, Inc.