

Stryker 1068 Eye Surgery Stretcher - Refurbished



- Dual Articulating Head Rest
- Collapsible Side Rails
- Central Brake and Steer
- \pm Height, \pm Trendelenburg
- Crank Operated Fowler
- 500 lbs Weight Capacity
- Includes New 3" Two Piece Pad
- Completely Refurbished

The Stryker 1068 Eye, Head, and Neck Surgical Stretcher features a dual articulating headrest for precision positioning of the patient during surgery.

The 1068 also has 8" casters for easy maneuverability, patient surface width of 26" and height range of 22" to 35.5" for easy/comfortable access for the surgeon.

The Stryker 1068 Eye, Head, and Neck stretcher comes completely refurbished with a comprehensive warranty and a new pad.

Specifications:

- Overall Length: 88"
- Overall Width: 31½"
- Patient Surface Width: 26"
- Height Range: 22.5" - 36"
- Weight Capacity: 500 lbs
- Fowler Angle: 0° - 90°
- Trend/Rev Trend: +17° / -15°
- Head Adjustment: \pm 3°
- Head Piece Tilt: \pm 45°
- Side Rails: 13" x 55"
- Wheel Diameter: 8"
- Pad Thickness: 3"