

Midmark 491 ENT Chair Refurbished



Stock Photo – Options May Vary

The Midmark 419 ENT/Otolaryngology chair makes exams and procedures easier with the ability to position patients from a low height of 22" to a maximum height of 31". The power backrest can be positioned from a fully upright 90° position to completely flat. The auto-return feature allows for a full return to home position so the patient can easily dismount.

This power chair has a foot control for operating movements, and a rotating base (330°). Also, it comes with new upholstery in your choice of colors.

The Midmark 491 comes completely refurbished with a comprehensive warranty and new upholstery.

- Modern Design
- Foot Control
- Electrically Operated
- +/- Height, +/- Backrest
- 22" Low Height
- Adjustable Headrest
- Auto-Return
- Armrests
- Rotating Base
- Completely Refurbished

Specifications:

- Minimum Height: 22" (55.9cm)
- Maximum Height: 31" (78.7cm)
- Width (Min): 24" (61cm)
- Width (Max): 26.5" (67.3cm)
- Rotation: 330°
- Electrical: Domestic: 115v, 12AMP

***Brochure and images are a general overview of the product.
Actual product options may vary.*