

Midmark 630 Power Exam Table

Refurbished



Stock Photo – Options May Vary

- Modern Design
- Hand Control
- Electrically Operated
- +/- Height, +/- Tilt
- Foot Up/Down and Back Up/Down
- Adjustable Headrest
- Hideaway Adjustable Stirrups
- Pull-Out Cysto Pan
- Table Paper Holder
- Table Paper Strap
- Accessory Side Rails
- Optional Foot Control
- Completely Refurbished

The Midmark 630 exam table/chair boasts a very low minimum height of 19". This makes it very easy for the patient to get on and off the chair safely and confidently. This power chair comes with a hand control for operating movements, an accessory side rail, a paper holder, and an articulating headrest.

The Midmark Ritter 630 comes with a pull-out drain pan at the end of the seat section and pull-out heel stirrups. The compact headrest provides better access to head and neck areas.

The Midmark Ritter 630 comes completely refurbished with a comprehensive warranty and new upholstery.

Specifications:

- Minimum Height: 19"
- Maximum Height: 40"
- Maximum Length:
 - Foot Extended: 72"
 - Foot & Headrest: 83"
- Width:
 - Seat: 28"
 - Back: 24" - 26"
- Weight Capacity: 450 lbs.
- Paper Roll Holder: 18" x 3.5"
- Electrical: 115 VACS

***Brochure and images are a general overview of the product.
Actual product options may vary.*