

Steris/Amsco Reliance 444 Washer/Disinfector - Refurbished



Stock Photo – Options May Vary

The Reliance 444 Single-Chamber Washer/Disinfector is a mechanical washer equipped with an Eagle 3000 Stage 3™ control system. The washer is designed with six (6) adjustable cycles:

- Instruments, Gentle, Utensils, Glassware, Plastic goods, Rubber goods
- Four (4) additional cycles are available for customized programming

Each pre-programmed cycle is equipped with: Pre-wash, Wash, Rinse, Thermal Rinse

The Amsco 444 comes fully refurbished with a comprehensive warranty.

Optional accessories and functions include: Manifold Rack, Loading/Transfer Cart, and Pure Rinse.

- Vertical Sliding Door
- Automatic Detergent Injection
- Tinted Tempered Glass Door
- Steam or Electric Operation
- Single or Double Door (Pass-Through)
- Dual Speed Pump (7.5 hp/ 2hp)
- Processes (8) full size instrument trays in 35 to 42 minutes
- Drying temperature range 150 degrees F (66° C) to 240°F (116° C)
- Double-walled insulated chamber
- Thermal Printer
- Up to 1.5 hrs continuous operation
- Includes Air Compressor (New)

Specifications:

- Exterior Dims: 42" x 80" x 32-5/8"
- Chamber Load Capacity: 26-1/2" x 26-1/4" x 26"
- Load Height: 31"
- Electrical Requirements: Electric Heated – 3 Phase - Call for Details
- Hot Water Inlet – 1/2" NPT, 5.5 GPM Min, 9.3 GPM Max, 15 dyn to 50 static PSI, 110° Max Temp
- Cold Water (If applicable) – 1/2" NPT, 10.3 GPM Min, 12.9 GPM Max, 30 dyn to 50 static PSI, 60°F Max Temp
- Recorder Paper: Thermal
- Vented Units – 3" ID Flexible Duct
- UL Listed

***Brochure and images are a general overview of the product. Actual product options may vary.*