

Steris/Amsco Flexmatic Dual Scrub Sink - Refurbished



Stock Photo – Options May Vary

The Steris Amsco Flexmatic dual station scrub sink is a heavy duty, hospital grade surgical scrub sink that has two independent bays. The water flow is controlled by knee actuation. The operator simply presses their knee against the stainless steel kick plate to start the water flow. Stopping the water flow is as easy as simply pressing the kick plate again.

The design is such that splashing is limited due to the stainless steel formed bottom of the sink. Rosebud spray heads create the perfect coverage of water and the gooseneck faucets direct the water in the ideal direction. Temperature is controlled independently at each station by a handle above the basin.

The Steris/Amsco Flexmatic scrub sink comes completely refurbished and with a comprehensive warranty.

- Stainless Steel Construction
- Heavy Duty
- Knee Operation (knee kick plate)
- 2 Station
- Rosebud Spray Head
- Formed Bottom Reduces Splashing
- No Electric Required
- Includes Z-Track and Mounting Brackets (2)
- Optional Divider Screen Available
- Completely Refurbished
- Comprehensive Warranty

Specifications:

- Dimensions: 64"W x 27.3"D x 39.5"H
- Height Off of Floor: 12" (recommended)
- Unit Weight: 280 lbs
- Cold Water Inlet: ½" NPT
- Hot Water – ½" NPT, 120-140° F
- Water Pressure Required: 20 to 50 psig

***Brochure and images are a general overview of the product. Actual product options may vary.*