

Steris / Amsco QDJ-04 Warming Cabinet - Refurbished



Stock Photo – Options May Vary

The Steris/Amsco QDJ 04 Warming Cabinet is a full size warmer that can heat blankets and/or liquids. This unit features separate top and bottom compartments which have independent temperature controls.

The QDJ-04 temperature range can be selected between 90° and 160°F. In case of overheating, the unit will alarm, both visually and audible, and shut the heating elements off.

The Steris/Amsco QDJ-04 comes completely refurbished with a comprehensive warranty and 2 shelves.

- Top and Bottom Heating
- Independent Controls
- Digital Displays
- Standard 120V Plug
- Stainless Steel Construction
- Over Temp Alarms
- Includes 2 Shelves
- Completely Refurbished

Specifications:

- Overall Dimensions:
30" W x 26.5" D x 74.75" H
- Interior Dimensions:
 - Upper
23.75" W x 22.87" D x 15.5" H
 - Lower
23.75" W x 22.87" D x 38.25" H
- Power: 120V 13.7Amps
- Approx Weight: 353 lbs
- Temp Range: 90° - 160°

***Brochure and images are a general overview of the product.
Actual product options may vary.*