

Stryker Gynnie 1061 OB-GYN Stretcher - Refurbished



Stock Photo – Options May Vary

The Stryker Gynnie 1061 OB-GYN Stretcher is perfect for emergency departments, labor and delivery, or gynecological procedures. The Stryker 1061 features a pneumatic assist fowler (backrest), folding leg section, and built-in stirrups.

Additional functions include Trendelenburg, reverse Trendelenburg, along with central brake and steer, glide-away side rails, and integrated O2 tank well.

The Stryker Gynnie 1061 OB-GYN stretcher comes completely refurbished with a comprehensive warranty and a new pad.

- Pneumatic Assist Fowler
- Glide-Away Side Rails
- Central Brake and Steer
- ± Height, ± Trendelenburg
- 500 lbs Weight Capacity
- Built-in O2 Tank Holder
- Includes New 4" Pad
- Built-in Stirrups
- Completely Refurbished

Specifications:

- Overall Length: 81½"
- Overall Width: 31½"
- Patient Surface Width: 26"
- Height Range: 22.5" – 36"
- Weight Capacity: 500 lbs
- Fowler Angle: 0° - 90°
- Trend/Rev Trend: ± 18°
- Side Rails: 13" x 55"
- Wheel Diameter: 8"
- Pad Thickness: 4"

***Brochure and images are a general overview of the product.
Actual product options may vary.*