

Stryker 1068 Eye Surgery Stretcher - Refurbished



Stock Photo – Options May Vary

- Dual Articulating Head Rest
- Collapsible Side Rails
- Central Brake and Steer
- \pm Height, \pm Trendelenburg
- Crank Operated Fowler
- 500 lbs Weight Capacity
- Includes New 3" Two Piece Pad
- Completely Refurbished

The Stryker 1068 Eye, Head, and Neck Surgical Stretcher features a dual articulating headrest for precision positioning of the patient during surgery.

The 1068 also has 8" casters for easy maneuverability, a patient surface width of 26", and a height range of 22" to 35.5" for easy/comfortable access for the surgeon.

The Stryker 1068 Eye, Head, and Neck stretcher comes completely refurbished with a comprehensive warranty and a new pad.

Specifications:

- Overall Length: 88"
- Overall Width: 31½"
- Patient Surface Width: 26"
- Height Range: 22.5" - 36"
- Weight Capacity: 500 lbs
- Fowler Angle: 0° - 90°
- Trend/Rev Trend: +17° / -15°
- Head Adjustment: \pm 3°
- Head Piece Tilt: \pm 45°
- Side Rails: 13" x 55"
- Wheel Diameter: 8"
- Pad Thickness: 3"

***Brochure and images are a general overview of the product.
Actual product options may vary.*