

Stryker 1069 Eye Surgery Stretcher - Refurbished



Stock Photo – Options May Vary

- Dual Articulating Head Rest
- Collapsible Side Rails
- Central Brake and Steer
- \pm Height, \pm Trendelenburg
- Gas Spring Assist Backrest
- 500 lbs Weight Capacity
- Includes New 3" Two Piece Pad
- Completely Refurbished

The Stryker 1069 Eye, Head, and Neck Surgical Stretcher offers the greatest head-end clearance in the industry. This stretcher features a dual articulating headrest for precision positioning of the patient during surgery and a gas spring-assisted backrest.

This stretcher is ideal for either the superior or temporal approach. Positioning controls tuck neatly within the headpiece and allow full articulation. The 1069 also features fold-down side rails, four IV receptacles, and four-wheel, steel-ring brakes with foot-end brake and steer pedals.

The Stryker 1069 comes completely refurbished with a comprehensive warranty and a new pad.

Specifications:

- Overall Length: 90"
- Overall Width: 31½"
- Patient Surface Width: 26"
- Height Range: 22.5" - 34"
- Weight Capacity: 500 lbs
- Fowler Angle: 0° - 90°
- Trend/Rev Trend: \pm 18°
- Head Adjustment: \pm 4.5°
- Head Piece Tilt: \pm 45°
- Side Rails: 13" x 55"
- Wheel Diameter: 8"
- Pad Thickness: 3"

***Brochure and images are a general overview of the product.
Actual product options may vary.*