

## GE Marquette MAC 1200 EKG Unit - Refurbished



Stock Photo – Options May Vary

The GE/Marquette Mac 1200 digital, 12-lead electrocardiograph system offers comprehensive ECG solutions with practical features that help meet the needs of hospitals, clinics, office based practices and clinical trials.

Its advanced algorithm capabilities, seamless connectivity to the CS and MUSE Cardiology Information Systems and easy-to-use features provide the quality and convenience users require.

Includes: Patient Cable, Clips, Pack of Paper, New Battery, Power Cord and Operators Manual.

The GE/Marquette Mac 1200 EKG comes completely refurbished with a comprehensive warranty.

- 12 Lead EKG
- Interpretation
- Full Page Printout
- LCD Display
- Preview Screen (Waveform Display)
- Full Keyboard
- Patient Annotation
- Memory (40 Resting ECGs)
- Portable and Lightweight
- Modes: Automatic, Manual, and Arrhythmia
- EMR Compatible

### Specifications:

- Display Type: LCD
- Display Resolution: 320 x 240
- Dims: 14.5"W x 12.6"D x 3.7"H
- Weight (unit only): 12.3 lbs
- Power Requirement: 120v
- Battery Charge Time: 4 Hours
- Battery Type: NiCad, 18v, 1.3Ah
- Printer Type: Thermal Array

*\*\*Brochure and images are a general overview of the product.  
Actual product options may vary.*