



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Biolight M70 Fingertip Oximeter- New



Stock Photo – Options May Vary

The Biolight M70 is an “all-in-one” integrated finger device with simple operation for immediate, real-time spot checks. The M70 is reliable with both accuracy and durability, can easily change direction of display, and it has an automatic power-off design.

Standard Parameters: Digital SpO₂, Pulse Rate, Pulse Bar, and Waveform.

Technical Specifications:

- Display: OLED
- Dimension: 2.3” (W) x 1.3 (H) x 1” (D)
- Weight: 0.13 lbs
- Power Supply: 2 AAA Alkaline Batteries

Biolight SpO₂ (Pulse Oximetry):

- Range: 0~100%
- Resolution: 1%
- Accuracy: ± 2%

Pulse Rate:

- Range: 25bpm -250bpm

Standard Accessory:

- 2 AAA Alkaline Batteries, Manual, Hang String

Biolight M700 Oximeter Patient Monitor- New



Stock Photo – Options May Vary

The Biolight M700 features: handhold design, compact, light for carrying and handling. The M700 also features bright, easy-to-read LED numerals, accurate measurement during low perfusion, anti-movement, and advanced digital processing technology.

Standard Parameters: Digital SpO₂, Pulse Rate, Pulse Bar, and Pleth Waveform.

Technical Specifications:

- Display: OLED, Bi-Color
- Dimension: 5½” (W) x 1” (H) x 2.9” (D)
- Weight: 0.27 lbs
- Power Supply: 4 AA Alkaline Batteries

Biolight SpO₂ (Pulse Oximetry):

- Range: 0~100%
- Resolution: 1%
- Accuracy: ± 2%

Pulse Rate:

- Range: 25bpm -250bpm

Standard Accessory:

- 4 AA Alkaline Batteries, Adult Reusable Sensor, Manual

Biolight M800 Oximeter Patient Monitor- New



Stock Photo – Options May Vary

The Biolight M800 features a 2.4” LED display that is compact and light for carrying and handling. The M800 functions include: automatic display rotation, two selectable working modes, spot check and continuous monitoring.

Standard Parameters: Digital SpO₂, Pulse Rate, Pulse Bar, and Pleth Waveform.

Technical Specifications:

- Display: 2.4” OLED, Bi-Color
- Dimension: 2.3” (W) x 4.8” (H) x 1” (D)
- Weight: 0.13 lbs
- Power Supply: 3 AA Alkaline Batteries

Biolight SpO₂ (Pulse Oximetry):

- Range: 0~100%
- Resolution: 1%
- Accuracy: ± 2%

Pulse Rate:

- Range: 25bpm -250bpm

Standard Accessory:

- 3 AA Alkaline Batteries, Manual,

***Brochure and images are a general overview of the product. Actual product options may vary.*

VENTUREMEDICAL.COM