



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Biolight M800 Handheld Pulse Oximeter - New



Stock Photo – Options May Vary

The Biolight M800 features a 2.4” LED display that is compact and light for carrying and handling. The M800 functions include automatic display rotation, two selectable working modes, spot check and continuous monitoring.

Standard Parameters: Digital SpO2, Pulse Rate, Pulse Bar, and Pleth Waveform.

Manufacturer Warranty: 2 Years, 90 Days for Accessories

Technical Specifications:

- Display: 2.4” OLED, Bi-Color
- Dimension: 2.3” (W) x 4.8” (H) x 1” (D)
- Weight: 0.13 lbs
- Power Supply: 3 AA Alkaline Batteries

Biolight SpO2 (Pulse Oximetry):

- Range: 0~100%
- Resolution: 1%
- Accuracy: $\pm 2\%$

Pulse Rate:

- Range: 25bpm -250bpm

Standard Accessory:

- 3 AA Alkaline Batteries, Manual, Reusable Finger Probe



Stock Photo – Options May Vary

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.