

## Blickman 7921D Warming Cabinet - Refurbished



*Stock Photo – Options May Vary*

The Blickman 7921D warming Cabinet is a cost effective solution for blanket and fluid warming needs in your medical facility, yet has modern features.

The 7921D warmer has digital controls for adjusting temperature. In case of overheating, the unit will alarm, both visually and audible, and shut the heating elements off.

The Blickman 7921D Warming Cabinet comes completely refurbished with a comprehensive warranty and 3 shelves.

- Single Compartment
- Digital Microprocessor
- Adjustable Digital Controls
- Digital Display
- Cost Effective
- Standard 120V Plug
- Stainless Steel Construction
- Over Temp Alarms
- Includes 3 Shelves
- Completely Refurbished

### Specifications:

- Overall Dimensions :  
24" D x 30" W x 72" H
- Internal Dimensions:  
22" D x 27" W x 50"H
- Power: 120V - 8.3 Amps
- Standard AC Wall Plug
- Temp Range: 90° - 160°
- UL Listed

*\*\*Brochure and images are a general overview of the product.  
Actual product options may vary.*