

## DateX Aestiva 5 Anesthesia Workstation - Refurbished



*Stock Photo – Options May Vary*

The Aestiva 5 Anesthesia Workstation offers an excellent combination of quality and performance at an affordable price. The 7900 ventilator offer a full range of operational parameters. The Aestiva 5 has been optimized for patient safety through patient centered engineering while also being ergonomically configured for the clinician's comfort.

The 7900 Smartvent offers pressure and volume control for superior ventilation. Low flow delivery techniques are allowed by the innovative compensation system, resulting in more consistent delivery of tidal volumes by adjusting automatically to variables by the machine and the patient.

The Aestiva 5 comes completely refurbished with a comprehensive warranty.

- Includes 7900 Smartvent Ventilator
- Integrated Absorber
- O<sub>2</sub> and Tidal Volume Monitoring
- Integrated Breathing Circuit
- Electronic PEEP
- One Motion from mechanical to Manual Mode
- Tidal Volume Compensation
- Low Flow Capable
- 2 Vaporizer Manifold
- Storage Drawers
- Lighted Strip
- Proven DateX-Ohmeda intuitive user interface, including anesthesia-dedicated menu logic and alarm philosophy, as in all System/5 monitors
- Compact Unit for places where space is a premium
- Built-in back-up battery to handle sudden power-down situations

### Specifications:

- Height (unit only): 53.4"
- Width (unit only): 36.6"
- Depth (unit only): 32.7"
- Weight (unit only): 340 Lbs
- Power Requirement: 120v

*\*\*Brochure and images are a general overview of the product.  
Actual product options may vary.*