

GE Dash 3000 Patient Monitor – Refurbished



Stock Photo – Options May Vary

The GE Dash 3000 Patient Monitor can easily be transported from room to room. The Monitor features compact size, lightweight, easy to use carry handle and built in battery, which make it portable and perfect in almost any environment.

The GE Dash 3000 offers a compact control panel and control knob, with its easy-to-use menu system that enables you to freeze, record, or perform other operations conveniently.

Included with the unit: ECG Cable (3 Lead), Reusable SpO2 finger sensor, NIBP Hose and Adult BP Cuff, Roll of Recorder Paper, New Batteries, Power Cord, and Users Manual.

The GE Dash 3000 patient monitor comes fully refurbished with a comprehensive warranty.

- Multi-Parameter Monitor
- Measures:
 - ECG, NIBP, SpO2, Temp
- 3 or 5 Lead ECG Standard
- Optional Onboard Thermal Printer
- Fully User Configurable Display
- 4 Auto-Set Alarm Functions
- 8.7" TFT Display
- Battery Operated w/ Built-In Charger
- Lightweight Design
- Carry Handle
- Mounting Options Available
- Full Range of Accessories Available
- Optional Invasive Pressure Available
- Optional Co2 Available
- Optional Agent Monitoring Available
- Completely Refurbished

Specifications:

- Dims: 10.24"H x 11.02" W x 7.87" D
- Weight: 12 lbs
- Display Size (Screen): 8.7"
- Battery: Rechargeable, SLA
- Charge Time: 4 Hrs Max
- Power Requirement: 120v

***Brochure and images are a general overview of the product. Actual product options may vary.*