

Maquet Alphastar 1132.01 Surgery Table - Refurbished



Stock Photo – Options May Vary

The Maquet Alphastar 1132 is an ultra-modern operating table and a great choice for almost any surgical setting. The 1132 features included: wired hand control, battery operation, automatic electric floor locks, and return-to-level function

Additional functions include: \pm back, \pm leg, \pm height, \pm tilt, \pm flex, trend/rev trend, and comes with North American side rails which fit most accessories.

The Maquet Alphastar 1132.01 comes completely refurbished with a comprehensive warranty and a new pad set.

- C-Arm Compatible Table Top
- Wired Hand Control
- 500 lbs Weight Capacity
- Automatic Electric Floor Locks
- Return-to-Level Function
- Battery Operated
- Accessory Side Rails
- Includes New Pad Set
- Completely Refurbished

Specifications:

- Overall Length: 82½"
- Overall Width: 22"
- Patient Surface Width: 20"
- Height Range: 27" - 44"
- Weight Capacity: 500 lbs
- Reverse Orientation: 300 lbs
- Trend/Rev Trend: \pm 30°
- Back Rest: + 68° / - 25°
- Leg Rest: + 80° / -105°
- Lateral Tilt: 18°
- Flex and Reflex Capable
- Power: 110-240 VAC
- Pad Thickness: 2"

***Brochure and images are a general overview of the product.
Actual product options may vary.*