



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Blickman 7922 Warming Cabinet - Refurbished



- Double Wall Insulated Construction
- Stainless Steel
- Magnetic Gasket Closure System
- Adjustable Temp
- Tempered Glass Door
- Standard 120V Plug
- Completely Refurbished

The Blickman 7922 single door warming cabinets or blanket warmers are 35-1/2" tall and suited for use under a counter-top or in limited spaces. May be used for heating and storage of solutions or as a blanket warmer. Built-in temperature data recorder. Warming Cabinet can sit on an existing cabinet, or on a Blickman mobile stand. Easy-to-use keypad with controls for setting internal temperature is located conveniently above door.

The Blickman 7922 comes completely refurbished with new hospital grade power cord and comprehensive warranty.

Specifications:

- Dims: 26.625" D x 30" W x 35.5" H
- Power: 120V 7.5 Amps
- Temp Range: 95° - 150°
- Interior Capacity: 9.04

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.