



www.venturemedical.com

Venture Medical ReQuip, Inc.  
6008 Bonacker Drive  
Tampa, FL 33610  
800-627-3215 | 813-740-8640  
Fax: 800-466-1251  
sales@venturemedical.com

## Conmed Excalibur Plus PC Electrosurgical Unit - Refurbished



- Digital Display
- Bipolar and Monopolar
- Multipurpose ESU
- High Power Output
- Accurate and Consistent Cutting
- Three Blend Modes
- ARM Grounding Fault Indicator
- Accurate and Consistent Cut Performance
- Low Voltage Coag Mode
- Monopolar Foot-Switch Included
- Bipolar Foot-Switch & Cart (Available)
- Completely Refurbished

The Conmed Excalibur Plus PC is the most powerful Isolated ESU on the market with 10,000 volts PTP in Spray COAG. It is the only ESU that allows the surgeon to adjust power settings from the surgical field with a standard ES pencil.

The Excalibur Plus PC Patented Controlled Output Voltage™ in CUT maintains minimum necessary voltage for smooth cutting through various impedance tissues with minimal thermal damage.

The Aspen Return Monitor (ARM™) is the most sensitive pad monitoring system available. It combines automatic feedback circuitry and visual displays to provide the most sensitive contact quality monitoring system available today. The ARM™ function continuously monitors the safety of the patient to ground System pad interface with a dual foil ground pad. It disarms the ESU and alerts the staff if a pad fault condition occurs.

Included with the unit: Monopolar Foot Switch, Power Cord, and Users Manual.

The Conmed Excalibur Plus PC comes completely refurbished with a comprehensive warranty.

### Specifications:

- Monopolar Cut: 300 Watts
- Monopolar Coag: 120 Watts
- Bipolar: 50 Watts
- Blend 1: 180 Watts
- Blend 2: 120 Watts
- Blend 3: 80 Watts
- Low Voltage Coag: 80 Watts
- Height (unit only): 6.25"
- Width (unit only): 15.9"
- Length (unit only): 22.75"
- Weight (unit only): 41 lbs
- Power Requirements: 120V, 8A