




www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

GE Marquette MAC 1200 EKG Unit - Refurbished



- 12 Lead EKG
- Interpretation
- Full Page Printout
- LCD Display
- Preview Screen (Waveform Display)
- Full Keyboard
- Patient Annotation
- Memory (40 Resting ECGs)
- Portable and Lightweight
- Modes: Automatic, Manual, and Arrhythmia
- EMR Compatible

The GE/Marquette Mac 1200 digital, 12-lead electrocardiograph system offers comprehensive ECG solutions with practical features that help meet the needs of hospitals, clinics, office based practices and clinical trials.

Its advanced algorithm capabilities, seamless connectivity to the CS and MUSE Cardiology Information Systems and easy-to-use features provide the quality and convenience users require.

Includes: Patient Cable, Clips, Pack of Paper, New Battery, Power Cord and Operators Manual.

The GE/Marquette Mac 1200 EKG comes completely refurbished with a comprehensive warranty.

Specifications:

- Display Type: LCD
- Display Resolution: 320 x 240
- Dims: 14.5"W x 12.6"D x 3.7"H
- Weight (unit only): 12.3 lbs
- Power Requirement: 120v
- Battery Charge Time: 4 Hours
- Battery Type: NiCad, 18v, 1.3Ah
- Printer Type: Thermal Array