



[www.venturemedical.com](http://www.venturemedical.com)

Venture Medical ReQuip, Inc.  
6008 Bonacker Drive  
Tampa, FL 33610  
800-627-3215 | 813-740-8640  
Fax: 800-466-1251  
[sales@venturemedical.com](mailto:sales@venturemedical.com)

## Midmark 540 Stretcher - Refurbished



- $\pm$  Height,  $\pm$  Trendelenburg
- Gas Spring Assist Backrest
- Fold Away Side Rails
- 4 Wheel Locking
- Central Brake and Steer
- Built-in O2 Tank Holder
- 350 lbs Weight Capacity
- Includes New 4" Pad
- Completely Refurbished

The Midmark 540 Stretcher is a versatile, full-functioning stretcher that is perfect for the operating room, emergency room, post-operative care unit or transport. The 540 has a 350 lbs weight capacity, central brake, four wheel locking and steering along with 8" casters for easy maneuverability and control.

Additional features include: pneumatic assist fowler 90°, hydraulic height adjustment (foot pedal on foot end of stretcher) 24" to 33", fold-down side rails and Trend/Reverse Trend.

The Midmark 540 Stretcher comes completely refurbished with a comprehensive warranty and a new pad.

### Specifications:

- Weight Capacity: 350 lbs.
- Overall Length: 83"
- Overall Width: 29.5"
- Height Range: 24" - 33"
- Fowler Backrest: 0° - 90°
- Trendelenburg: +12°
- Reverse Trendelenburg: -7°
- Patient Surface: 77" x 27"
- Side Rails: 15" x 63-5/8"
- Casters: 8"
- Pad Thickness: 4"