



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Skytron 3500 Slide-Top Operating Table - Refurbished



- C-Arm Compatible Table Top
- Wired Hand Control
- Return-to-Level Function
- Full Slide-Top
- 500 lbs Weight Capacity
- 180° Table Top Rotation
- Electro-Hydraulic Floor Locks
- North American Side Rails
- Includes New 3 Piece Pad

The Skytron 3500 Slide-Top surgical table is a multi-purpose operating table that features the ability to perform full body imaging. This table is a great choice for vascular, neuro, urology and orthopedics.

The 3500 slide-top features include: a wired hand control, 21" top slide capability, electro-hydraulic break, C-Arm compatibility, and return-to-level.

Additional functions include: \pm back, \pm leg, \pm height, \pm tilt, trend/rev trend, and comes with North American side rails which fit most accessories.

The Skytron 3500 slide-top comes completely refurbished with a comprehensive warranty and a new pad set.

Specifications:

- Overall Length: 79"
- Overall Width: 21.75"
- Patient Surface Width: 19.75"
- Height Range: 26" - 45"
- Weight Capacity: 500 lbs
- Trend/Rev Trend: + 25° / -20 °
- Back Rest: + 90° / - 40°
- Leg Rest: + 2° / -98°
- Lateral Tilt: \pm 20°
- Slide-Top: 21"
- Power: 120 V, 300 Watts
- Pad Thickness: 2"