

## Steris/Amsco 3085 SP Surgery Table - Refurbished



The Steris/Amsco 3085SP operating table is a great choice for most surgical procedures. The table has all standard functions  $\pm$  back,  $\pm$  leg,  $\pm$  height,  $\pm$  tilt,  $\pm$  flex, trend/rev trend. The 3085SP also features a radiolucent tabletop, return-to-level function which is actuated from the wired hand control, is battery operated, and has a 1000 lbs weight capacity in normal position, and 500 lbs in reverse orientation.

The wide range of accessories available is an added benefit. Accessories available include: arm boards, stirrups, anesthesia screens, shoulder positioners, urology extensions, orthopedic leg holders and extensions, and much more.

The Steris/Amsco 3085SP surgery table comes completely refurbished with a comprehensive warranty and includes a new 3 piece pad set.

- Wired Hand Control
- Accessory Side Rails
- Radiolucent Top
- Electro-Hydraulic Brakes
- $\pm$  Back,  $\pm$  Leg,  $\pm$  Height,  $\pm$  Trendelenburg,  $\pm$  Tilt,  $\pm$  Flex
- Integrated Kidney Lift
- Return-to-Level Function
- Battery Operated
- 1,000 lbs Weight Capacity (Normal Position)
- Includes New 3 Piece Pad
- Completely Refurbished

### Specifications:

- Overall Length: 76"
- Overall Width: 21  $\frac{5}{8}$ "
- Patient Surface Width: 20"
- Height Range: 27" - 44"
- Weight Capacity: 1000 lbs
- Trend/Rev Trend:  $\pm$  25°
- Back Rest: + 55° / - 25°
- Leg Rest: + 80° / -105°
- Side Tilt:  $\pm$  18°
- Flex: 20° Back  $\downarrow$  / 25° Seat  $\downarrow$
- Reflex: 25° Back  $\uparrow$  / 35° Seat  $\uparrow$
- Side Rails: 13" x 55"
- Power: 120 Volts, 4.5 Amps
- Pad Thickness: 2"