



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Steris Amsco Quantum Surgical Lights -Refurbished



- Flexible, Versatile Lighting
- Excellent Depth of Illumination (26")
- Automatic Re-Lamping
- Color temperature 4200k
- Includes 360 Degree Swivel Arm Support
- Intensity: 10,000 ft candles (107,600 lux)
- Wall Mounted Intensity Controller
- Color Rendering Index: 92
- Focusable Light Head
- Excellent Shadow Control
- Quartz Halogen Bulbs
- Refurbished

The Amsco Quantum surgical lighting system is designed to provide exceptional illumination in almost any surgical setting. The 24" light head is easily maneuvered and the auto re-lamping feature means that you should never be without light, even if a bulb fails during surgery.

The large, pre-focused beam of light is pre-focused, and yields an incredible 10,000 foot candles (107,600 lux). Impressive natural color rendition at 4200 Kelvin and a color rendering index (CRI) of 92 provides are also key features of this system.

The Steris Amsco Quantum surgical light can come in a single or dual head configuration. Also included is the wall controller and transformer.

The surgical lights come fully refurbished and includes a warranty.

Specifications:

- Halogen Bulb: 22 V / 220 W
- Intensity Pattern (large): 7.5"
- Intensity Pattern (small): 6.75"
- Clearance (rotation): 14.5'
- Intensity: Up to 10,000 Ft Candles
- Diameter of Light Head: 24"
- Power Requirement (volts/Hz):
 - 110-120V, 50/60 Hz
 - 230-240V, 50/60 Hz
- Power Requirements (amps):
 - 7A (dual head)
 - 3.5A (single head)

VENTUREMEDICAL.COM

© 2007-2025 Venture Medical ReQuip, Inc.