



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Steris/Amsco QDJ-03 Warming Cabinet - Refurbished



- Tabletop Model
- Space Saving Design
- Digital Display
- Standard 120V Plug
- Stainless Steel Construction
- Over Temp Alarms
- Includes 1 Shelf
- Completely Refurbished

The Steris/Amsco QDJ-03 series of tabletop blanket warmers is a perfect option for your department or surgery center. The QDJ-03 is a cost-saving solution versus the cost of a comparable new unit.

Features Include: Single Compartment, Stainless Steel Construction, and Quiet Operation, Overheat Alarm Digital Display.

The Steris/Amsco QDJ-03 comes with a rack, completely refurbished and with a comprehensive warranty.

Specifications:

- Overall Dimensions:
30" W x 26½" D x 24.25" H
- Inner Dimensions:
24" W x 23" D x 15½" H
- Power: 120V 6.8Amps
- Approx Weight: 260 lbs
- Temp Range: 90° - 160°
- Voltage: 110v