

## Ritter/Midmark M9 Ultraclave Autoclave - Refurbished



*Stock Photo – Options May Vary*

- Automatic/Electronic Autoclave
- Overheat Protection
- Visual and Audible Alarms
- Automatic Door Opener for Venting
- Includes Trays
- 9” Diameter Chamber
- Optional Printer
- Completely Refurbished
- Comprehensive Warranty

The Ritter/Midmark M9 Ultraclave Autoclave is one of the best-selling tabletop automatic/electronic autoclaves in the industry. The M9 features a 9” diameter chamber with automatic venting.

The Midmark/Ritter M9 offers easy, effective sterilization. It couldn’t be easier to use because it lets you load it, set it and forget it. Select your cycle and the time is set for you. Press “start” and the M9 fills with the correct amount of water. Once the process is complete, the door opens automatically to dissipate steam and give you the fastest drying time possible.

The Midmark/Ritter M9 Ultraclave Autoclave comes completely refurbished with a comprehensive warranty and includes trays.

### Specifications:

- Overall Length: 18.81”
- Overall Width: 15”
- Countertop Area: 18” Depth
- Usable Volume: 3.5 Gallons
- Trays: Two Large, Two Small
  - 12” L x 7.31” W
  - 12” L x 5.62” W
- Chamber Size: 9” x 15” Depth
- Electrical Requirements:
  - 120 Volt, 12 Amps
- Recorder Paper: 2.25” W
- Weight: 70 lbs.
- UL and CUL Listed

*\*\*Brochure and images are a general overview of the product.  
Actual product options may vary.*