


## NorthStar Beach Chair POS-2000 - Refurbished



*Stock Photo – Options May Vary*

The Northstar POS-2000 Beach Chair positioner for shoulder procedures provides an easy to use and safe method for shoulder procedures. The unit includes a tiltable head rest and drops panel design to expose the operative side of patient. The back rest is adjustable by pulling a trigger latch on back of unit. Adjusts easily with gas spring assist. The Northstar POS 2000 has a 450 lbs weight capacity and mounts to almost any surgery table with accessory side rails.

Included with the unit: Base Unit, Tilting Headrest, Lateral Support and Pad, Side Rails (removable), New Pads and Straps.

The NorthStar Beach Chair comes fully refurbished with a comprehensive warranty.

- Safe and Convenient
- Easy to Use
- Universal Mounting Brackets
- Fits Almost Any Surgery Table
- Gas Spring Assist Back Rest
- Adjustable Side Support
- Tilting Head Rest
- Drop Panel Design
- 450 Lbs Capacity
- New Pads and Straps
- Completely Refurbished

### Specifications:

- Overall Dims: 39" x 19.75"
- Side Panel: 26.5" x 6.5" (ea)
- Weight: 38 lbs
- Load Capacity: 450 Lbs
- Patient Height – No Restrictions
- Articulation Range: 0° - 90°

*\*\*Brochure and images are a general overview of the product. Actual product options may vary.*