

Physio-Control Lifepak 12 Biphasic Defibrillator With AED and Pacing - Refurbished



Stock Photo – Options May Vary

The Physio-Control Lifepak 12 is known for its quality and reliability. The Lifepak 12 comes equipped with industry standard monitoring tools for adult and pediatric patients.

The Lifepak 12 features ECG, manual defibrillation, noninvasive pacing, AED functions, and biphasic technology.

Included with the unit: Multi-Function Cable, ECG Cable and Leadwires, New Batteries (2), Operator's Manual, 1 Set of Adult Pads and Recorder Paper.

Optional Accessories Available: AC Adapter (Charges Directly from Wall Outlet, No Separate Battery Charger Needed!) and External Paddles.

The Physio-Control Lifepak 12 comes completely refurbished with a comprehensive warranty.

- Manual Defibrillation
- Pacing
- Large Screen
- Durable Design
- 3 Lead ECG
- AED Function
- Biphasic Technology
- Hands-Free Multi-Function Cable
- Recorder
- Alarms
- Code Summary
- Lightweight and Compact
- Completely Refurbished

Specifications

- Max Joules: 360j
- Charge Time: Less than 10 Seconds
- Display Size (Screen): 5.5"W x 4.2"H
- Battery: Dual Battery Capability
- Recorder Paper: 50mm
- Height: 12.5"
- Width: 15.6"
- Depth: 9.1"
- Weight: 14.8 lbs
- Power Requirement: 90-132VAC,47-63Hz

***Brochure and images are a general overview of the product.
Actual product options may vary.*