

Skytron 6500HD Hercules Bariatric Surgical Table - Refurbished



Stock Photo – Options May Vary

The Skytron 6500HD operating table is a modern multi-purpose surgical table and a great choice for almost any surgical setting. The 6500HD features include wired hand control, 180° tabletop rotation, hydro-electric brakes, C-Arm compatibility, and 850 lbs. weight capacity in normal position.

The Skytron 6500HD has all standard functions ± back, ± leg, ± height, ± tilt, ± flex, trend/rev trend, and comes with North American side rails which fit most accessories.

This Skytron 6500HD comes completely refurbished with a comprehensive warranty and a new pad set.

- C-Arm Compatible Table Top
- Wired Hand Control
- 850 lbs Weight Capacity
- 180° Table Top Rotation
- Electro-Hydraulic Floor Locks
- Stainless Steel Base and Column
- North American Side Rails
- Swivel Casters
- Includes New 3 Piece Pad

Specifications:

- Overall Length: 76"
- Overall Width: 21 ¾"
- Patient Surface Width: 20"
- Height Range: 26" - 47"
- Weight Capacity: 850 lbs
- Trend/Rev Trend: ± 30°
- Back Rest: + 90° / - 40°
- Leg Rest: + 20° / -105°
- Lateral Tilt: 30°
- Flex and Reflex Capable
- Power: 120 Volts
- Pad Thickness: 2"

***Brochure and images are a general overview of the product. Actual product options may vary.*