

Steris/Amsco 3085 SP Surgery Table - Refurbished



Stock Photo – Options May Vary

The Steris/Amsco 3085SP operating table is a great choice for most surgical procedures. The table has all standard functions \pm back, \pm leg, \pm height, \pm tilt, \pm flex, trend/rev trend. The 3085SP also features a radiolucent tabletop, return-to-level function which is actuated from the wired hand control, is battery operated, and has a 1000 lbs weight capacity in normal position, and 500 lbs in reverse orientation.

The wide range of accessories available is an added benefit. Accessories available include: arm boards, stirrups, anesthesia screens, shoulder positioners, urology extensions, orthopedic leg holders and extensions, and much more.

The Steris/Amsco 3085SP surgery table comes completely refurbished with a comprehensive warranty and includes a new 3 piece pad set.

- Wired Hand Control
- Accessory Side Rails
- Radiolucent Top
- Electro-Hydraulic Brakes
- \pm Back, \pm Leg, \pm Height, \pm Trendelenburg, \pm Tilt, \pm Flex
- Integrated Kidney Lift
- Return-to-Level Function
- Battery Operated
- 1,000 lbs Weight Capacity (Normal Position)
- Includes New 3 Piece Pad
- Completely Refurbished

Specifications:

- Overall Length: 76"
- Overall Width: 21 $\frac{5}{8}$ "
- Patient Surface Width: 20"
- Height Range: 27" - 44"
- Weight Capacity: 1000 lbs
- Trend/Rev Trend: \pm 25°
- Back Rest: + 55° / - 25°
- Leg Rest: + 80° / -105°
- Side Tilt: \pm 18°
- Flex: 20° Back \downarrow / 25° Seat \downarrow
- Reflex: 25° Back \uparrow / 35° Seat \uparrow
- Side Rails: 13" x 55"
- Power: 120 Volts, 4.5 Amps
- Pad Thickness: 2"

***Brochure and images are a general overview of the product. Actual product options may vary.*