

Steris/Amsco Century V120 Sterilizer - Refurbished



Stock Photo – Options May Vary

The Steris Amsco Century V120 sterilizer offered by Venture Medical is a Prevac unit, which is designed for fast, efficient sterilization of surgical instruments and related items. This sterilizer boasts a large 20" x 20" chamber size with a depth of 38". The unit is user-friendly and has enhanced ease of use with its intuitive touch screen interface. The vertical sliding door, activated by a foot pedal, makes loading and unloading the unit much easier than with older models that used a traditional swinging door.

All refurbished units from Venture Medical come standard with a new 30kW steam generator, unless house steam is available. The Steris Amsco V120 is available in freestanding or recessed mounting configurations and includes one loading rack and two shelves.

The Steris Amsco V120 sterilizer comes fully refurbished with a comprehensive warranty.

- Vertical Sliding Door
- 20" x 20" x 38" Chamber Size
- Touch Screen Display
- Loading Rack and Two Shelves
- Prevac Operation
- Multiple Preprogrammed Cycles
- Self Tests and Alarms
- Help Screen for Troubleshooting
- Steam or Electric Operation
- Single or Double Door (Pass-Through)
- Freestanding or Recessed Mounting
- Optional Transfer/Loading Cart
- Thermal Printer

Specifications:

- Exterior Dims: 45"D x 30"W x 74.5"H
- Chamber Size: 20"x 20" x 38"D
- Control Panel Electrical Requirements: 120 Volt, 1-Phase, 9.5 Amps
- Steam Generator Electrical Requirements: Electric Heated – 3 Phase - Call for Details
- Water Feed: 1" NPT, 30-50 psi
- Drain: 1-1/2" ODT
- Operating Temperature: 250° - 270°F
- Recorder Paper: 2.25" W
- UL Listed

***Brochure and images are a general overview of the product. Actual product options may vary.*